



Speedboat™
Notch

Multi-modal endoscopic device with advanced tissue traction features

Advanced Energy-powered bipolar RF dissection and SHF microwave coagulation for gastric, oesophageal, duodenal and colonic procedures in the Upper and Lower GI tracts



Anything is Possible
with the Right Approach



CREO
MEDICAL

Speedboat Notch: Safety & Stability

Advanced Tissue Traction

Speedboat Notch is engineered with advanced tissue traction, precise tissue dissection, and an innovative notch tip design for depth perception —advancements that provide precise control, improved tissue interaction, and optimal procedural outcomes.

Seamlessly integrating with CROMA advanced bipolar RF cutting technology and super high-frequency (SHF) 5.8GHz microwave coagulation, Speedboat Notch offers a versatile, all-in-one solution for complex tissue resection and efficient haemostasis, for Lower GI, Upper GI Gastric, and Upper GI Esophageal procedures.

Speedboat Notch Features:


- Advanced tissue traction provides precision and control
- Precise tissue dissection provides controlled cutting capability
- Innovative notch tip design for depth perception, engineered to provide safe dissection and coagulation for haemostasis





Lower GI




Upper GI Gastric




Upper GI Esophageal



Innovative Notch Tip Design for Depth Perception



Precise Tissue Dissection



Advanced Tissue Traction

Speedboat Notch: Precision & Control


Precise advanced bipolar RF cut

The blade design controls the depth of penetration and provides a focused pathway of energy delivery at lower voltage <460 V. The voltage/current is adjusted based on tissue impedance automatically to maintain power density for a smooth, high quality and precise cut.

Adaptive waveform automatically adjusts parameters to tissue and **balances coagulation** during cutting to minimise bleeding..

Designed to perform safe, precise, contact cut with clean margins to provide **high quality histology** samples and promote healing.

Protective Hull protects the muscle bed from unwanted thermal injury by maintaining a constant distance from the energy source, allowing cutting close to the muscle bed.




On-demand, controlled SHF microwave coagulation

SHF 5.8 GHz microwave energy distributes heat evenly across the treatment area. The microwave energy is delivered through the tip of the instrument, allowing for application of energy through the distal end.

5.8 GHz Microwave enables **controlled penetration depth**, unaffected by tissue resistance, **reducing the risk of perforation and charring¹**.

Fast delivery of coagulation with controlled spread and depth of penetration to prevent or treat bleeding immediately without changing devices.



CROMA Advanced Energy Platform: The power behind Speedboat Notch

The CROMA Advanced Energy Platform precisely controls **advanced bipolar RF** and **SHF 5.8 GHz microwave** energy to enable a suite of flexible endoscopic devices designed to deliver:

- A unique **usability and safety profile¹⁻⁶**
- Optimal **tissue effect¹⁻⁶**
- Improved **clinical and economic outcomes⁶**
- **Expanded capabilities** in therapeutic endoscopy

Find out more:





Specifications

Specification	Speedboat Notch (Long)	Speedboat Notch (Short)
Product Reference	PRD-SB1-003	PRD-SB1-004
Min. Channel Size		2.8mm
Max Catheter Size		2.4mm
Working Length / Full Length	1.9m / 2.3m	1.25m / 2.3m
Advanced Bipolar RF (Cut)		15 -35 Watts
Super High Frequency Microwave (Coag)		08 - 10 Watts

Visit: www.creomedical.com for more information

References

1. Data on file
2. Microwave coagulation of blood vessels during advanced colonoscopic polypectomy: first results in humans. Zacharias P. Tsiamoulos et al. published in United European Gastroenterology Journal; 2016; 2 (Supplement 1). [https://www.giejournal.org/article/S0016-5107\(17\)31361-5/pdf](https://www.giejournal.org/article/S0016-5107(17)31361-5/pdf)
3. A new approach to endoscopic submucosal tunneling dissection: the "Speedboat-RS2" device. Zacharias P. Tsiamoulos et al. published in Endoscopy. <https://www.thieme-connect.de/products/ejournals/html/10.1055/a-0875-3352>
4. Endoscopic submucosal tunneling dissection: use of a novel bipolar radiofrequency and microwave-powered device for colorectal endoscopic submucosal dissection. Thomas R. McCarty, Hiroyuki Aihara. Published in Video GIE, official video journal of the American Society of Gastrointestinal Endoscopy. [https://www.videogie.org/article/S2468-4481\(20\)30090-4/fulltext](https://www.videogie.org/article/S2468-4481(20)30090-4/fulltext)
5. Tsiamoulos et al. First results using Speedboat Tunneling technique in colorectal submucosal dissection – clinical outcomes and procedure time prediction models. Poster presented at UEG 2020. <https://ueg.eu/library/first-results-using-speedboat-tunneling-technique-in-colorectal-submucosal-dissection-clinical-outcomesandprocedure-time-prediction-models/240928>
6. Cost-effectiveness analysis of Speedboat submucosal dissection in the management of large non-pedunculated colorectal polyps, based on 50 patients. Authors: Amir Ansari-pour, Mehdi Javanbakht, Adam Reynolds, Zacharias Tsiamoulos. Data on file.

Creo Medical Ltd.

Unit 2, Creo House
Beaufort Park Way
Chepstow
NP16 5UH, UK
+44 (0) 1291 637 300
customerservice@creomedical.com

Creo Medical Inc.

100 Reserve Road
Suite B400
Danbury,
CT 06810, USA
+1 866-226-1170
UScustomerservice@creomedical.com

Creo Medical Pte Ltd.

Creo Medical Pte Ltd.
8 Commonwealth Lane
#04-03C
Singapore
149555
Customer Service
+1 866-226-1170