

# Smart Bulb

Simple and intuitive



# Smart Bulb

Smart Bulb is designed to improve visibility by efficiently clearing away debris in the field-of-view during ureterorenoscopy.

## In Line Flow Switch (IFS)

The In Line Flow Switch (IFS) allows comfortable, single-handed on/off control of flow and prevents inadvertent spillage of the fluid.

## Dual Flow

- A continuous, low intensity background flow regulated by adjusting the height of the irrigation fluid reservoir, and;
- An accelerated flow to clear debris or move the position of the stone when the bulb is compressed.

## Tubing Extension

- A tubing extension has been added to each unit. This can be attached to the Smart Bulb to make working with a flexible scope easier.
- Allows for the scrub nurse to assist in the procedure comfortably.

PRODUCT CODE	DESCRIPTION	QTY
URO 037	Smart Bulb	10



Wismed CC  
71 Princess Alice Avenue,  
Durban, South Africa, 4001  
Tel: +27 (0) 31 206 0177  
orders@wismed.co.za



TNMC Devices Limited  
Ivana Lucica 2a/19  
10000 Zagreb, Croatia  
Tel: +38 59 1622 9658

