


# EchoHood

## Ballonnet pour Écho-Endoscope

**EchoHood** is a balloon designed to withstand the demands of echo-endoscopy while guaranteeing a very high quality ultrasound image.

**A complete range,**  
constantly following the  
evolution of echo-endoscopes.







STERILE







# EchoHood

## Balloons for Echo-Endoscopes

### EchoHood – Balloons for Olympus™ Echo-Endoscopes





Pictogram	Reference	Description	Opening(s)	QTY
	AS1070000	Balloon for echo-duodenoscopes	2 openings	20
	AS1460000	Balloon for echo-endoscopes	2 opening	20
	AS1830000	Balloon for echo-duodenoscopes	1 opening	20
	AS2230000	Balloon for echo-endoscopes	2 openings	20
	AS3280000	Balloon for echo-bronchoscopes	2 openings	20
	AS3680000	Balloon for echo-duodenoscopes	1 opening	20

### EchoHood – Balloons for Pentax™ Echo-Endoscopes

Pictogram	Reference	Description	Opening(s)	QTY
	AS1470000	Balloon for echo-duodenoscopes	1 opening	20
	AS2070000	Balloon for echo-duodenoscopes	2 openings	20
	AS3710000	Balloon for echo-duodenoscopes	2 openings	20
	AS3940000	Balloon for echo-bronchoscopes	2 openings	20

AS1470000 is supplied with fixing rings and a mounting plug.

### EchoHood – Balloons for Fujifilm™ Echo-Endoscopes

Pictogram	Reference	Description	Opening(s)	QTY
	AS2660000	Balloon for echo-endoscopes	2 openings	20
	AS2800000	Balloon for echo-endoscopes	1 opening	20
	AS3560000	Balloon for echo-endoscopes	1 opening	20
	AS3610000	Balloon for echo-endoscopes	2 openings	20

## Characteristics

- ▶ **Constant thickness membrane** for excellent image quality.
- ▶ **Rigid and opaque individual packaging** to protect the balloon from U.V. and shocks.
- ▶ **Individual sterile case** to avoid any risk of contamination by contact.
- ▶ **Balloons individually checked** for quality assurance.
- ▶ **Correspondence table “ Echo-endoscopes - Balloons”** available on request.



**Description:** Echo-Hood Balloon for echo-endoscopes. **Indications:** Digestive and bronchial echo-endoscopy. **Class:** Is. **Notifying body:** CE0459. **Manufacturer:** Prince Medical. **Proper use:** Before use, carefully read the instructions in the package insert or on the label of the device. **Version:** EH\_E\_2. Rev.02 Date 08-06-23

Creo Medical authorised local distributor

**Creo Medical UK Ltd**

Chestnut House, N.Minster Business Pk.

YO26 6QR York, United Kingdom

**Email:** info.uk@creomedical.com

**Tel:** +44 (0) 1904 786888

Technical changes and mistakes reserved. Images may show non standard options.

