

# SmartScale - Aros - SmartFlow

## Wireless flowmeters

- ▶ Non Invasive Diagnostics.
- ▶ Several solutions for different needs.



## SmartScale - PC Flow

### The wireless Flowmeter with PC and network integration

#### Advantages of the SmartScale

- ▶ Wireless – easy to use and clean
- ▶ Fully automatic.
- ▶ Excellent measurement quality.
- ▶ Fixed or variable scale.
- ▶ Automatic analysis including Liverpool, Siroky and paediatric nomograms.
- ▶ Easy to use software with SQL database.
- ▶ High security level for data protection.
- ▶ Hygienic - designed to prevent overflow and splashes.
- ▶ Robust and durable.



98800002

SmartScale



## Aros - Flow-EMG

### Combined EMG and Flow system

#### Advantages of Aros

- ▶ Aros has all the features and advantages of the SmartScale.
- ▶ Wireless synchronisation with the EMG.
- ▶ Rechargeable EMG connected to the computer USB
- ▶ Optional Biofeedback software with skin electrodes, anal and vaginal probes.
- ▶ Possibility to combine with our residual volume measurement system, Pinit.

98400000

Aros

## SmartFlow - Printer Flow

### The "classic approach" with wireless design and unique extras

#### Advantages of the SmartFlow:

- ▶ Wireless, portable and easy to use.
- ▶ Touch screen control.
- ▶ High measurement quality especially in small flows
- ▶ Fully automatic, power-safe mode, Auto-zero, Auto-start, Auto-stop, Auto-print.
- ▶ Fast printing (80mm thermal paper).
- ▶ Automatic analysis and measurements.
- ▶ Data security, stores the latest 4 studies within internal memory to allow re-printing.
- ▶ Hygienic - designed to prevent overflow and splashes.
- ▶ Optional: wall mounted.
- ▶ Optional: PC link allows studies to be saved into an SQL database allowing integration with the hospital information system.



00004440

SmartFlow

**Description:** Urine Flowmeters. **Indications:** To perform Mictiometry studies with and without EMG signal. **Class:** IIa. **Notifying Body:** SGS. **Manufacturer:** Creo Medical S.L. **Caution:** Please read carefully the Instruction for Use Manual prior to start using the equipment. **Version:** FM\_E\_3 Rev.03 Date 23-05-23

Creo Medical authorised local distributor

#### Creo Medical UK Ltd

Chestnut House, N.Minster Business Pk.

YO26 6QR York, United Kingdom

Email: info.uk@creomedical.com

Tel: +44 (0) 1904 786888

Technical changes and mistakes reserved. Images may show non standard options.

