

ScopeNet

Net Snare



Safe recovery
of polyps,
foreign bodies
and food
impactions.



CREO
MEDICAL


ScopeNet

Net Snare

Reference	Ø Catheter (mm)	Catheter length (cm)	Net size I x L (mm)	QTY
FST-FG-18-160-C2	1,8	160	25 x 55	10
FST-FG-23-230-C3	2,3	230	30 x 60	10
FST-FG-23-230-C4	2,3	230	40 x 75	10

Characteristics

- ▶ **360° rotatable handle** for easy positioning.
- ▶ **Ergonomic design of the scissor-type handle** that facilitates one-handed operation.
- ▶ **Teflon sheath** that provides **excellent glide** in the operating channel.
- ▶ **Expandable net** for the recovery of large foreign bodies.
- ▶ **Preserves specimen morphology** during extraction.
- ▶ **Sterile packaging** ensures a product free of contamination.



Description: SCOPENET net extraction device for polyps and foreign bodies. **Indications:** This product is used with an endoscope to remove objects, food boluses or polyps from the human body after resection with a polypectomy loop. **Class:** IIa. **Notified body:** CE0197. **Manufacturer:** Foresight Technology Co, Ltd. **Proper use:** Before use, carefully read the instructions in the package insert or on the device label. **Version:** NS_E_2 Rev.02 Date 15-07-2022

Creo Medical authorised local distributor

Creo Medical UK Ltd

Chestnut House, N.Minster Business Pk.

YO26 6QR York, United Kingdom

Email: info.uk@creomedical.com

Tel: +44 (0) 1904 786888

Technical changes and mistakes reserved. Images may show non standard options.

