

ScopeWire


Guide Wire

A wide range of guidewires for selective bile duct cannulation, placement and exchange of ERCP devices.

- ▶ **5 cm hydrophilic end.**
- ▶ **Two-tone zebra coating** allows correct visualization of correctly visualise the movement of the guide wire on the endoscopic image.
- ▶ **2 diameters: 0.025" and 0.035".**
- ▶ **2 lengths: 260 cm and 450 cm.**
- ▶ **Standard or rigid model.**
- ▶ **Straight or angled** hydrophilic end.

**A wide range of
guidewires for
selective bile duct
cannulation, placement
and exchange of ERCP
devices.**

ScopeWire - Guide Wire					
Référence	Ø (Inch / mm)	Length (cm)	Rigidity	Type flexible end	QTY
GW-25265041	0,025 / 0,63	260	Standard	Straight	10
GW-25455041	0,025 / 0,63	450	Standard	Straight	10
GW-35265041	0,035 / 0,89	260	Standard	Straight	10
GW-35455041	0,035 / 0,89	450	Standard	Straight	10
GW-25265043	0,025 / 0,63	260	Rigid	Straight	10
GW-25455043	0,025 / 0,63	450	Rigid	Straight	10
GW-35265043	0,035 / 0,89	260	Rigid	Straight	10
GW-35455043	0,035 / 0,89	450	Rigid	Straight	10
GW-35455042	0,035 / 0,89	450	Standard	Angled	10



Description: ScopeWire guide. **Indications:** Placement of devices in the biliary and digestive tracts by endoscopy. **Class:** IIa. **Notified body:** CE0123. **Manufacturer:** Hangzhou AGS Medtech Co, Ltd. **Proper use:** Before use, carefully read the instructions in the package insert or on the device label. **Version:** SGUIDE_E_5. Rev.05. Date 15-07-2022