


UroCatch


Extraction Basket

4 references where the extraction basket is declined:

- ▶ In 2 lengths: 90cm and 120cm.
- ▶ In 3 diameters: 1.9Fr, 2.2Fr and 3.0Fr.
- ▶ With or without tip
- ▶ Straight or helical wire basket.



To perform endoscopic removal of upper urinary tract stones.




UroCatch

Extraction Baskets

Reference	Ø Catheter (Fr)	Length (cm)	Basket	Ø Panier (mm)	Tip	QTY
USB-30090-001	3,0	90	4 helicoidal wires	13 - 16	✓	1
USB-30090-009	3,0	90	4 straight wires	13 - 16	✓	1
USB-19120-042	1,9	120	4 straight wires	11 - 16	✗	1
USB-22120-041	2,2	120	4 straight wires	11 - 16	✗	1

Characteristics

- ▶ Ergonomic handle, easy to **remove** and **re-assemble**.
- ▶ **360° rotatable** basket.
- ▶ **Radiopaque** sheath and basket.
- ▶ **Nitinol basket**: shape memory material, resistant to deformation.
- ▶ Models **without Tip** to **facilitate the capture of stones**, especially in the calcareous funds.



UC_F_1

Description: UroCatch stone extraction basket. **Indications:** Indicated for choledocholithiasis and kidney stones. **Class:** IIa. **Notified body:** CE1984. **Manufacturer:** Plasti-Med Plastik Medikal Ürünler San.ve Tic.Ltd.şti. **Proper use:** Before use, carefully read the instructions in the package insert or device labelling. **Version :** UC_F_1. Rev.01 Date 22-12-2021

Distributed in UK by:

Creo Medical UK Ltd.

Chestnut House, Northminster Business Park
YO26 6QR York, United Kingdom

Email: info.uk@creomedical.com

Tél: +44 (0) 1904 786888

