


ScopeClip^P

Refillable haemostatic clip

1 pre-mounted clip + 2 or 4 additional clips



**High-performance
Economical & Eco-
friendly**




Rechargeable



360° rotating



Repositionable



**CREO
MEDICAL**

ScopeClip^P

Refillable Hemostatic Clip

Reference	Number of clips	Clip opening (mm)	Clip angle (°)	Ø Catheter (mm)	Catheter length (mm)	Colour	CDT
HC-2311-2-13555504	1 + 2	11	135	2,6	2300	Green	10
HC-2316-2-13555504	1 + 2	16	135	2,6	2300	Blue	10
HC-2311-4-13555504	1 + 4	11	135	2,6	2300	Green	10
HC-2316-4-13555504	1 + 4	16	135	2,6	2300	Blue	10

Also available with a 90° clip angle, please contact us.

Features

- ▶ **The only refillable clip that is both rotatable and repositionable**, both the pre-assembled clip and the refills:
 - **Rotates 1:1 360° in both directions** for easy positioning on the target area.
 - **Repositionable**: can be reopened as often as necessary to correct positioning.

- ▶ **Cost-effective** compared with conventional single-use clips.

- ▶ **Quick and easy to recharge.**

Reloading takes no more time than opening a new clip! See for yourself how efficient reloading is by flashing this QR Code :


Scan here




[Click here](#)

- ▶ **And above all, ScopeClip^P is eco-friendly:**

ScopeClip³
1 unit = 3 units of conventional clips



ScopeClip⁵
1 unit = 5 units of conventional clips



Less than :

- ▶ **Raw materials**
- ▶ **Packaging and paper**: 1 instruction leaflet/box
- ▶ **Transport**: for the same quantity of clips
- ▶ **Waste**: divided by 3 or 5 for the same number of clips fitted
- ▶ **Space-saving**: simplified stock management

Name: ScopeClip^P single-use refillable haemostasis clip. **Indications:** Endoscopic insertion of clips for haemostasis or marking. **Class:** IIb. **Notified body :** CE0123. **Manufacturer:** Hangzhou AGS MedTech Co. Hangzhou AGS MedTech Co, Ltd. **Correct use:** Before use, carefully read the instructions for use or labelling of the device. **Reimbursement:** Product not included in the LPPR. **Version:** SCL3_F_2. Rev.02. Date 25-03-24

Distributed in France by:

Creo Medical SASU

9, avenue Jean Prouvé

88101 Saint-Dié des

Vosges

Email info.fr@creomedical.com

Tel: + 33 3 29 58 45 45 / 46 46

