


ScopeBite

Digestive Biopsy Forceps







A wide range of biopsy forceps designed to **optimize tissue sampling**, even in hard-to-reach areas.



ScopeBite

Single-use digestive biopsy forceps



ScopeBite - Biopsy forceps standard model

	Reference	Catheter Ø (mm)	Length (cm)	Jaws	Spike	Jaw opening (mm)	Withdrawal volume (mm ³)	Color coding	PKG
	GBF-18165011	1.8	160	Smooth	✗	4.5	≥ 3.0	Grey	10
	GBF-23165011	2.3	160	Smooth	✗	6.7	≥ 7.3	Orange	10
	GBF-23235011	2.3	230	Smooth	✗	6.7	≥ 7.3	Blue	10
	GBF-23165012	2.3	160	Smooth	✓	6.7	≥ 7.3	Orange	10
	GBF-23235012	2.3	230	Smooth	✓	6.7	≥ 7.3	Blue	10
	GBF-23165013	2.3	160	Toothed	✗	6.7	≥ 7.3	Orange	10
	GBF-23235013	2.3	230	Toothed	✗	6.7	≥ 7.3	Blue	10
	GBF-23165014	2.3	160	Toothed	✓	6.7	≥ 7.3	Orange	10
	GBF-23235014	2.3	230	Toothed	✓	6.7	≥ 7.3	Blue	10

ScopeBite - Biopsy forceps with swinging jaws

	Reference	Catheter Ø (mm)	Length (cm)	Jaws	Spike	Jaw opening (mm)	Withdrawal volume (mm ³)	Color coding	PKG
	GBF-23165021	2.3	160	Smooth	✗	6.7	≥ 7.3	Orange	10
	GBF-23235021	2.3	230	Smooth	✗	6.7	≥ 7.3	Blue	10
	GBF-23165022	2.3	160	Smooth	✓	6.7	≥ 7.3	Orange	10
	GBF-23235022	2.3	230	Smooth	✓	6.7	≥ 7.3	Blue	10

ScopeBite - Biopsy forceps large capacity with swinging jaws

	Reference	Catheter Ø (mm)	Length (cm)	Jaws	Spike	Jaw opening (mm)	Withdrawal volume (mm ³)	Large capacity	Color coding	PKG
	GBF-23165023	2.4	160	Toothed	✗	9.0	≥ 9.6	✓	Orange	10
	GBF-23235023	2.4	230	Toothed	✗	9.0	≥ 9.6	✓	Blue	10
	GBF-23165024	2.4	160	Toothed	✓	9.0	≥ 9.6	✓	Orange	10
	GBF-23235024	2.4	230	Toothed	✓	9.0	≥ 9.6	✓	Blue	10

Characteristics

- ▶ **2 models available:**
 - Standard model,
 - Large capacity swinging jaws model.
- ▶ **Two types of jaws:**
 - Smooth with sharp edges,
 - Toothed.
- ▶ **Large capacity model with swinging jaws:**
 - Larger **sampling volume ≥ 9.6 mm³**,
 - Double fenestration of each jaw,
 - **Allows the jaws to pivot** at the site to be biopsied for accurate tissue sampling, even in hard-to-reach areas.

Denomination: Single-use biopsy forceps. **Intended use:** Sampling of living tissue from the digestive tract and respiratory tract under endoscopy. **Class:** IIa. **Notified body:** CE0344. **Manufacturer:** Transmed (China) Co., Ltd. **Proper use:** Before use, carefully read the instructions on the instructions leaflet or device labeling. **Version:** SBMT_E_2. Rev.02. Date 07-08-2025

Creo Medical UK Ltd

Creo House, Unit 2 Beaufort Park Way,
Chepstow, NP16 5UH, UK
Email: info.uk@creomedical-europe.com
Tel: +44 (0) 1904 786888

Creo Medical Inc

100 Reserve Road, Suite B400
Danbury, CT 06810, USA
Email: uscustomerservice@creomedical.com
Tel: +1 866-226-1170

Creo Medical Pte Ltd

8 Commonwealth Lane, #04-03C
Singapore, 149555
Email: info.apac@creomedical.com

