

ScopeClip 360


Single-use Hemoclip



Repositionable



Rotatable



Enhanced bidirectional rotation ensures pinpoint accuracy when positioning on the lesion.




CREO
MEDICAL

ScopeClip 360

Hemostatic clip

ScopeClip 360 - Hemostatic clip **rotatable** and **repositionable**

Reference	Clip span (mm)	Clip angle	Catheter Ø (mm)	Working length (mm)	Color code	PKG
HC-23111355104	11	115°	2.6	2300	Green	10
HC-231613551044	16	115°	2.6	2300	Blue	10



Characteristics

- ▶ **Enhanced bidirectional rotation**, ensuring accurate and controlled positioning on the target area.
- ▶ **Repositionable**: can be reopened to adjust positioning.
- ▶ **Once opened, the clip can be closed to a smaller opening width** adapting to the lesion or the anatomy of the target area if its maximum opening is too wide.
- ▶ The clip's opening and closing mechanism keeps it securely closed without continuous handle activation.



We also have a single-use reloadable, rotatable and repositionable clip.

For more information on this model, please contact us.

Denomination: Single-use hemostatic clip. **Intended Use:** This device is intended to be used for hemostasis for mucosal/submucosal defects in digestive tract. **Class:** IIb. **Notified body:** CE0344. **Manufacturer:** Transmed (China) Co., Ltd. **Proper use:** Before use, carefully read the instructions on the instruction leaflet or device labeling. **Version:** SCMT_E_1. Rev.01. Date 29-04-2025

Creo Medical UK Ltd

Creo House, Unit 2 Beaufort Park Way,
Chepstow, NP16 5UH, UK

Email: info.uk@creomedical-europe.com

Tel: +44 (0) 1904 786888

Creo Medical Inc

100 Reserve Road, Suite B400
Danbury, CT 06810, USA

Email: uscustomerservice@creomedical.com

Tel: +1 866-226-1170

Creo Medical Pte Ltd

8 Commonwealth Lane, #04-03C
Singapore, 149555

Email: info.apac@creomedical.com

