

ScopeInject

Injection Needle

6 configurations where the needle is available in :

- 3 diameters: 21, 23, and 25 G,
- 2 lengths: 4 and 6 mm.




ScopInject Injection Needle

Reference	Catheter		Needle		PKG
	Ø (mm)	Length (cm)	Ø (G / mm)	Length (mm)	
IN-241621045081	2.3	160	21 / 0.8	4.0	10
IN-242321045081	2.3	230	21 / 0.8	4.0	10
IN-241621065081	2.3	160	21 / 0.8	6.0	10
IN-242321065081	2.3	230	21 / 0.8	6.0	10
IN-241623045082	2.3	160	23 / 0.6	4.0	10
IN-242323045082	2.3	230	23 / 0.6	4.0	10
IN-241623065082	2.3	160	23 / 0.6	6.0	10
IN-242323065082	2.3	230	23 / 0.6	6.0	10
IN-241625045083	2.3	160	25 / 0.5	4.0	10
IN-242325045083	2.3	230	25 / 0.5	4.0	10
IN-241625065083	2.3	160	25 / 0.5	6.0	10
IN-242325065083	2.3	230	25 / 0.5	6.0	10

Characteristics

- ▶ **Triple bevelled needle** for **safe and precise punctures**.
- ▶ Transparent sheath for visualizing any needle movement.



Denomination: Injection needle. **Intended Use:** Endoscopic submucosal injection in digestive tract and injection in the treatment of esophageal and gastric varices. **Class:** IIa. **Notified body :** CE0344. **Manufacturer:** Transmed (China) Co., Ltd. **Proper use:** Before use, carefully read the instructions on the instruction leaflet or device labeling. **Version :** SIMT_E_1. Rev.01. Date 30-04-25

Creo Medical UK Ltd

Creo House, Unit 2 Beaufort Park Way,
Chepstow, NP16 5UH, UK
Email: info.uk@creomedical-europe.com
Tel: +44 (0) 1904 786888

Creo Medical Inc

100 Reserve Road, Suite B400
Danbury, CT 06810, USA
Email: uscustomerservice@creomedical.com
Tel: +1 866-226-1170

Creo Medical Pte Ltd

8 Commonwealth Lane, #04-03C
Singapore, 149555
Email: info.apac@creomedical.com

