

ScopeSnare COLD

Single Use Cold Snare



A thin 0.23 mm
diameter cutting wire
for a clean sectioning of
the lesion.



ScopeSnare COLD


Single Use Cold Snare

ScopeSnare COLD - Rotatable Cold Polypectomy Snare

	Reference	Ø Catheter (mm)	Length (cm)	Snare Ø (mm)	Snare Shape	PKG
	CPS-2410235336	2.3	230	10	Diamond	10
	CPS-2415235336	2.3	230	15	Diamond	10
	CPS-2410235335	2.3	230	10	Oval	10
	CPS-2415235335	2.3	230	15	Oval	10

Characteristics

- ▶ **Rotatable**, to facilitate positioning of the snare on the target lesion.
- ▶ **Multi-strand braided cutting wire with a diameter of 0.23 mm**, allowing a clean resection of the lesion.
- ▶ **“Diamond” snare shape**, specifically engineered to maximize opening width for accurate placement.
- ▶ **Blue handle coding** indicates compatibility with cold resection techniques.
- ▶ Snare opening width marked on the handle for ease of use.



Braided multi-strand Ø 0.23 mm



Denomination: Single Use Cold Snare. **Intended Use:** Endoscopic resection of diminutive polyps, sessile polyps, pedunculated polyps and tissue from within the gastrointestinal tract. **Class:** IIa. **Notified body:** CE0344. **Manufacturer:** Transmed (China) Co., Ltd. **Proper use:** Before use, carefully read the instructions on the instruction leaflet or device labeling. **Version:** SSCMT_S_1 Rev.01. Date 29-04-2025

Creo Medical UK Ltd

Creo House, Unit 2 Beaufort Park Way,
Chepstow, NP16 5UH, UK
Email: info.uk@creomedical-europe.com
Tel: +44 (0) 1904 786888

Creo Medical Inc

100 Reserve Road, Suite B400
Danbury, CT 06810, USA
Email: uscustomerservice@creomedical.com
Tel: +1 866-226-1170

Creo Medical Pte Ltd

8 Commonwealth Lane, #04-03C
Singapore, 149555
Email: info.apac@creomedical.com

